



STABILITY MAPPING USING THE TSI FORMULATION OPTIMIZATION

Context

Suspensions are unstable colloidal systems that undergo many destabilization phenomena (sedimentation, flocculation, creaming, etc.) that can be due to various causes (viscosity, particle size, attractive forces, etc.). Therefore, it may be complex for the formulator to determine the origin of these processes to be able to overcome them and to get a stable product.

One of the strategies to avoid these destabilization phenomena is to increase the viscosity of the formulation. The purpose of this paper is to determine the effect of increasing polymer concentration in a simple dispersion of TiO₂ in different mixtures of glycerol and water.

Measurement protocol

16 samples were studied, adjusting the formulation with various ratios:

- Ratio glycerol/water: 0 / 0.33 / 0.66 / 1
- Xanthan concentration: 0 / 0.05 / 0.1 / 0.2
- 0.50% wt. TiO₂

Identification of phenomena

Thanks to the TURBISCAN device, destabilization phenomena were identified in all formulations.




Figure 1. Delta Backscattering; water/glycerol: 2/1; xanthan: 0

The graph above can be separated into 3 parts to understand the phenomena:



- At the top (3), the BS signal decreases because the concentration decreases;
- In the middle (2), the BS signal evolves, because of the size variation;
- At the bottom (1), the BS signal increases because the concentration increases;

Qualification of phenomena

It is possible to rank the stability of the samples, depending on the composition ([xanthan], [glycerol/water]), thanks to the TSI as shown below.

TSI sums all the variations detected in the samples in terms of size and/or concentration. The higher is the TSI, the worse is the stability. As shown below, it is possible to plot a stability diagram, representing the stability (TSI), as a function of the product formulation.




Figure 2. Effect of Composition on TSI



As expected, the stability of these dispersions is improved by:

- Increasing the viscosity of the continuous phase, simply by adjusting the ratio between glycerol and water;
- Adding xanthan to the continuous phase, is another way to increase its viscosity, thus slowing down the migration of the particles.

It is also noticeable that optimal stability is obtained for the samples containing more than 0.1% of xanthan and more than 66% of glycerol in the continuous phase.

Conclusion

The TURBISCAN enables to choose the best formulation in terms of stability. This is possible in every kind of concentrated dispersion without any dilution. The TURBISCAN is also a unique tool for monitoring the kinetics of destabilization of the system and quantifying it precisely after identifying the phenomena.



TURBISCAN RANGE

The world leader in Stability & Dispersibility Analysis