

ACCURATE ANALYSIS OF PORE STRUCTURES FOR ACTIVATED CARBON

PORE SIZE DISTRIBUTION WITH A NEW GCMC KERNEL BASED ON A SLIT-PORE MODEL WITH CARBON SURFACE HETEROGENEITY

SUITABLE APPLICATIONS



Catalyst



Adsorbents



Battery cell




Electronic component

APPLICABLE MATERIALS



Carbon




New carbon slit GCMC kernel with roughness surface

MEASUREMENT & ANALYSIS


- Materials:** Activated carbon fiber (FT300-25; Kurarey Chemical) pretreated at 305°C for 3h in vacuum
- Ads. isotherm:** In N₂ at 77.4K, Ar at 87.3K by BELSORP MAX series (Microtrac)
- PSD Analysis:** GCMC (smooth), GCMC (roughness), QSDFT

RESULTS

ADSORPTION ISOTHERM FITTING



PORE SIZE DISTRIBUTION



I GCMC with smooth model

An ideal adsorption isotherm with a smooth kernel exhibits partially different behavior from the experimental data, resulting in a valley in the pore size distribution at porediameters of approximately 0.9 nm.

I GCMC with roughness model **NEW**

The ideal adsorption isotherm with the updated kernel with surface roughness best fits the experimental isotherm and represents a more accurate pore size distribution.

I QSDFT

The PSDs derived from the QSDFT method also showed systematic valleys and peaks, indicating insufficient results.

EQUIPMENT



Gas & Vapor Adsorption Analyzers

BELSORP MAX X SERIES

Measurement items

- BET surface area
- Pore size distribution
- Gas & vapor adsorption



High-precision gas adsorption measurement system

BELSORP MAX G

Measurement items

- BET surface area
- Pore size distribution
- Gas adsorption



The new **GCMC kernel with roughness surface** and **QSDFT kernel** are implemented in the analysis software BELMASTER (Ver. 7.4.3.1~) and are suitable for the data from **BELSORP MAX X SERIES** and **BELSORP MAX G**.

MICROTRAC

PARTICLE CHARACTERIZATION

Microtrac Inc.
11 Penns Trail • Newtown, PA 18940 • USA
Phone: +1 866 473 8724 • marketing@microtrac.com

Microtrac Retsch GmbH
Retsch-Allee 1-5 • 42781 Haan • Germany
Phone: +49 2104 2333 300 • info@microtrac.com

MicrotracBEL Corp.
8-2-52 Nanko Higashi, Suminoe-ku • Osaka 559-0031 • Japan
Phone: +81 6-6655-0362 • international@microtrac-bel.com

Microtrac Formulaction SAS
3-5 rue Paule Raymondis • 31200 Toulouse • France
Phone: +33 (0)5 62 89 29 29 • contact.fr@mtf.verder.com

www.microtrac.com

part of **VERDER**