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



Automatic optimization of conditions  
such as gas introduction volume

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#### WATER BATH

Easily swap the Dewar with a water bath for measurement temperatures ranging from  $-10^{\circ}\text{C}$  to  $70^{\circ}\text{C}$ . A refrigerated / heated circulator is required.

#### HEATER

In-situ pretreatment of the sample from  $50^{\circ}\text{C}$  to  $450^{\circ}\text{C}$ . With this option, samples can be pretreated directly at the measurement port, without the need to transfer the sample cell before analysis. Particularly useful for sensitive samples.

#### GAS SELECTOR

The gas selector enables the connection of up to 4 adsorptives simultaneously (1 x Helium + 4 x adsorptive).

HIGH-PRECISION GAS ADSORPTION ANALYZER BELSORP MAX G

## BELCONTROL: NEW OPERATION SOFTWARE

The versatility of BELSORP equipment is truly world leading. The numerous features and capabilities are complemented by BELCONTROL the intuitive and user-friendly operation software. It guides the user step-by-step through the analysis process. This includes the setup of analysis conditions, executing the measurements, when to fill and setup the liquid nitrogen or other bath, when to replace the gas cylinder, the degassing steps, and much more. The software is designed to make the instrument accessible and operable to everyone, including inexperienced users.

For inexperienced users or for measurements of unknown samples, BELCONTROL only requires basic sample information (name, mass, etc.), pre-treatment conditions (if not performed externally) and the measurement range.

Detailed control of the configuration and measurement settings is possible to optimize the measurement conditions (e.g. dosing settings, equilibrium criteria, leak test option, etc.). This allows the user to fully customize the sample analysis to his needs.

| BELCONTROL              |             |                                                                           |
|-------------------------|-------------|---------------------------------------------------------------------------|
| Quick BET               | tak         | Multi-point BET surface area in less than 20 minutes                      |
| Helium-Free Measurement | tak         | AFSM <sup>TM</sup> 2 enables He-free measurements with unmatched accuracy |
| Adsorption Kinetics     | opcjonalnie | Rate of adsorption measurements for diffusion analysis                    |

### MORE BELCONTROL FEATURES

- | Overlaying adsorption / desorption isotherms and comparing the measured data between the various ports during measurement
- | All pressures, temperatures, valve actuations, etc. are stored in trend data, allowing for immediate examination
- | A system check function is available for

diagnosing the instrument status

- | E-mail notification automatically transmits measurement status and results
- | Interactive program in Japanese or English ensures easy, reliable operation
- | Extensive help functions, including step-by-step instructions during operation


HIGH-PRECISION GAS ADSORPTION ANALYZER BELSORP MAX G

## **OPERATION SOFTWARE**

The versatility of the BELSORP instruments, with their many features and possibilities, is complemented by our intuitive and user-friendly operation software. It guides the user step-by-step through the analysis process including setting up the analysis conditions, executing the measurements, when to fill and setup the liquid nitrogen or other bath, when to replace the gas cylinder, the degassing steps, and much more. The software has been crafted to make the instrument accessible to everyone, even inexperienced users.



Two measurement modes are offered: 'Simple Mode' and 'Professional Mode'. Simple Mode enables easy operation, requiring minimal setup. You only need to enter the basic sample information (name, mass, etc.) and then select your pretreatment conditions and analysis conditions from a list of preset options. This is ideal for inexperienced users or measurements of unknown samples. If a prior measurement with comparable sorption behavior has been performed, GDO can be used to reduce the measurement time. Professional Mode allows detailed control of the configuration and measurement settings (e.g., dosing settings, equilibrium criteria, leak-check option, etc.), allowing the user to completely customize sample analysis.


GAS ADSORPTION MEASUREMENTS OF POROUS & NON-POROUS MATERIALS: MOFS, ZEOLITES, CARBONS,  
AND MORE

## **EVALUATION SOFTWARE BELMASTER**

A sorption isotherm is defined as the relationship between the adsorbed amount onto an adsorbent and the equilibrium pressure of the adsorptive gas – commonly related to the saturation vapor pressure – at constant temperature. The gas sorption isotherm (e.g. nitrogen) delivers information about the specific surface area, pore size distribution and pore volume of measured material. In the following graph some exemplary sorption isotherms are shown.


The specific surface area (SSA) refers to the accessible surface area of the sample and is of great importance in adsorption, homogeneous and heterogeneous catalytic reactions. The specific surface area can be calculated according to ISO 9277 with the BET method (BET: Brunauer, Emmett and Teller) or Langmuir method. The following graph shows an exemplary evaluation of the specific surface area according to the BET method in our BELMASTER software:

Select the right pressure range (multi point BET) or a measurement point (single point BET) and the surface area will be calculated automatically. Further, our BELMASTER software also provides the BET surface area calculation according to ISO 9277 Annex C (also known as Rouquerol-plot), which is recommended for microporous materials.

The outstanding feature of the BELSORP MAX G is the availability of three different transducers (1000, 10 and 1 / 0.1 torr), enabling very low pressure adsorption isotherms, starting from  $p/p_0 = 10^{-8}$  ( $N_2$  @ 77K) up to atmospheric pressure. As a result, pore size distributions from 0.35 up to 500 nm can be obtained. Further, low surface area measurements down to  $0.0005 \text{ m}^2/\text{g}$  by using krypton gas are possible.


The characterization of low specific surface area materials such as non-porous metallic materials, glass substrates, and films can be difficult or impossible with traditional gases like nitrogen (77 K) and argon (77 K or 87 K) due to detection limits. Instead, krypton gas adsorption can be used at liquid nitrogen temperature to determine the BET specific surface area down to  $0.0005 \text{ m}^2/\text{g}$ .

## HIGH-PRECISION GAS ADSORPTION ANALYZER BELSORP MAX G

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
## HIGH-PRECISION GAS ADSORPTION ANALYZER BELSORP MAX G

**ZASADA DZIAŁANIA****Analysis of Pore Size Distribution via NLDFT & GCMC Method**

The classical pore size distributions (PSD) are calculated by INNES method (slit shape) and BJH, DH, CI methods (cylinder shape), which evaluate mesopores based on the capillary condensation theory. HK (slit), SF (cylinder), and CY (cage) methods can also be used to evaluate micropores based on the adsorption potential theory. The DA method, as well as the DR method, are also commonly used for pore volume evaluation and pore structure evaluation. The new PSD and capacity evaluation methods, NLDFT and GCMC (optional for BELSORP MINI X), are used for wide range analyses (from micro- to meso and macropores) and more accurate methods, specified in ISO15901-2.

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[www.microtrac.com/belsorp-max-g](http://www.microtrac.com/belsorp-max-g)